


# CUSTOMFLEX® ARTIFICIALIRIS

Case and Images provided courtesy of Kamran M. Riaz, MD, Dean McGee Eye Institute



## Patient's History:

- 69 y.o. patient with history of rock injury to left eye at age 6 with traumatic cataract removal, was left aphakic and with significant iris loss 360 degrees (functional aniridia)
- 2016: PPV with scleral fixated IOL (Gore-Tex® (off label use)) OS; CF visual acuity
- 2017: underwent PKP OS for corneal decompensation; CF visual acuity
- 2019: repeat PKP for failure of first PKP
- 2020: Though VA has been stable at CF for years, pt is a farmer and very much bothered by glare/photophobia from OS that limits his outdoor time.

## Implantation/Post-Op Results:

- 2020: Underwent artificial iris implant surgery OS with CUSTOMFLEX® ARTIFICIALIRIS (with fiber)
- Though vision remains CF, pt is able to work in fields again with minimal sunglasses

