



## CUSTOMFLEX® ARTIFICIALIRIS: Photo of the Month

This patient is a 43 year old male with a past ocular history of congenital aniridia, glaucoma and visually significant cataracts with severe zonulopathy. His pre operative BCVA was 20/400 OU, and he could no longer work or see his two children clearly. His visual potential was excellent without significant nystagmus or foveal hypoplasia. His final vision after cataract extraction, IOL implantation and CUSTOMFLEX® ARTIFICIALIRIS was 20/40 and 20/30, with medically controlled glaucoma.



To submit your case, please contact us at:

513-872-1330  
schristensen@veo-ophthalmics.com

Images and case provided courtesy of Nicole Fram, MD, Advanced Vision Care, Los Angeles, CA.